

NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Quantity Available
E2611L	Gibson Assembly® Master Mx	50 reactions	2
E2611S	Gibson Assembly® Master Mx	10 reactions	2
E2621S	NEBuilder® HiFi DNA Assembly Master Mx	10 reactions	2
E2621X	NEBuilder® HiFi DNA Assembly Master Mx	250 reactions	2
E2621L	NEBuilder® HiFi DNA Assembly Master Mx	50 reactions	2
NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Quantity Available
M0202L	T4 DNA Ligase	100,000 units	2
M0202S	T4 DNA Ligase	20,000 units	2
M0202T	T4 DNA Ligase	20,000 units	2
M0202M	T4 DNA Ligase	100,000 units	2
NEB #	Markers & Ladders	Size	Quantity Available
N3200L	1 kb Plus DNA Ladder	1,000 gel lanes (0.5 ml)	2
N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	2
NEB #	Nucleic Acid Purification	Size	Quantity Available
T1010L	Monarch® Plasmid Miniprep Kit	250 preps	2
NEB #	PCR, Polymerases & Amplification Technologies	Size	Quantity Available
M0484L	OneTaq® Hot Start 2X Master Mix with Standard Buffer	500 reactions (50 µl)	2
M0484S	OneTaq® Hot Start 2X Master Mix with Standard Buffer	100 reactions (50 µl)	2
M0492S	Q5® High-Fidelity 2X Master Mix	100 reactions (50 µl)	2
M0492L	Q5® High-Fidelity 2X Master Mix	500 reactions (50 µl)	3
M0494L	Q5® Hot Start High-Fidelity 2X Master Mix	500 reactions	2
M0314L	RNase Inhibitor, Murine	15,000 units	2
M0314S	RNase Inhibitor, Murine	3,000 units	2
NEB #	RNA Synthesis and Modification	Size	Quantity Available
M0239L	T4 RNA Ligase 2 (dsRNA Ligase)	750 units	2