

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0367S	Blunt/TA Ligase Master Mix	50 reactions	M0367S
M0370S	Instant Sticky-end Ligase Master Mix	50 reactions	M0370S
M2200L	Quick Ligation™ Kit	150 reactions	M2200L
M2200S	Quick Ligation™ Kit	30 reactions	M2200S
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0202L	T4 DNA Ligase	100,000 units	M0202L
M0202M	T4 DNA Ligase	100,000 units	M0202M
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	N3200S
N3200L	1 kb Plus DNA Ladder	1,000 gel lanes (0.5 ml)	N3200L
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
N3012S	Lambda DNA HindIII Digest	150 gel lanes (0.3 ml)	N3012S
N3026S	ΦX174 DNA-HaeIII Digest	50 gel lanes (0.05 ml)	N3026S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
N0446S	Deoxynucleotide (dNTP) Solution Set	25 μmol	N0446S
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	M0210S
M0368S	ProtoScript® II Reverse Transcriptase	4,000 units	M0368S
M0491S	Q5® High-Fidelity DNA Polymerase	100 units	M0491S
M0297S	RNase H	250 units	M0297S
M0203S	T4 DNA Polymerase	150 units	M0203S
M0254S	Vent® DNA Polymerase	200 units	M0254S
NEB #	Restriction Endonucleases	Size	Stockroom #
R0161S	AccI	1,000 units	R0161S
R0652S	AfeI	200 units	R0652S
R0136S	BamHI	10,000 units	R0136S
R3136S	BamHI-HF®	10,000 units	R3136S
R0143S	BglI	2,000 units	R0143S
R0144L	BglII	10,000 units	R0144L
R0144S	BglII	2,000 units	R0144S
R0197S	Clal	1,000 units	R0197S
R0101S	EcoRI	10,000 units	R0101S

R3101S	EcoRI-HF [®]	10,000 units	R3101S
R0195S	EcoRV	4,000 units	R0195S
R0103S	HincII	1,000 units	R0103S
R0104S	HindIII	10,000 units	R0104S
R3104S	HindIII-HF [®]	10,000 units	R3104S
R3142S	KpnI-HF [®]	4,000 units	R3142S
R0191S	NarI	500 units	R0191S
R0193S	NcoI	1,000 units	R0193S
R0193L	NcoI	5,000 units	R0193L
R0189S	NotI	500 units	R0189S
R3189S	NotI-HF [®]	500 units	R3189S
R0140S	PstI	10,000 units	R0140S
R3140S	PstI-HF [®]	10,000 units	R3140S
R0151S	PvuII	5,000 units	R0151S
R3156S	SacI-HF [®]	2,000 units	R3156S
R0157S	SacII	2,000 units	R0157S
R0138S	SalI	2,000 units	R0138S
R3138S	SalI-HF [®]	2,000 units	R3138S
R0141S	SmaI	2,000 units	R0141S
R0182S	SphI	500 units	R0182S
R0187S	StuI	1,000 units	R0187S
R0145L	XbaI	15,000 units	R0145L
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0180S	XmaI	500 units	R0180S
R0194S	XmnI	1,000 units	R0194S
NEB #	RNA Synthesis and Modification	Size	Stockroom #
S1404S	m ⁷ G(5')ppp(5')G RNA Cap Structure Analog	1 μmol	S1404S

NEB Representative: David Gdula, Ph.D.
gdula@neb.com
(978) 380-2970

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com