

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0226S	CpG Methyltransferase (M.SssI)	100 units	M0226S
M2200S	Quick Ligation™ Kit	30 reactions	NE2200S
M0202S	T4 DNA Ligase	20,000 units	NE202S
M0202T	T4 DNA Ligase	20,000 units	NE202CS
M0392S	β-Agarase I	100 units	NE392S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	NE323-2S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	NE323-1S
N3233S	Low Molecular Weight DNA Ladder	100 gel lanes (0.1 ml)	NE323-3S
NEB #	Nucleic Acid Purification	Size	Stockroom #
T1120S	Monarch® Spin DNA Gel Extraction Kit	50 preps	T1120S
T1110L	Monarch® Spin Plasmid Miniprep Kit	250 preps	T1110L
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 μmol	N0447S
N0446S	Deoxynucleotide (dNTP) Solution Set	25 μmol	NE446
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	NE210S
M0203S	T4 DNA Polymerase	150 units	NE203S
M0267S	Taq DNA Polymerase with ThermoPol® Buffer	400 units	NE267S
M0254S	Vent® DNA Polymerase	200 units	NE254S
NEB #	Restriction Endonucleases	Size	Stockroom #
R0117S	AatII	500 units	R0117S
R3552S	Agel-HF®	300 units	R3552S
R0114S	Apal	5,000 units	NE114S
R0558S	Ascl	500 units	NE558S
R0152S	Aval	2,000 units	NE152S
R0174S	AvrII	100 units	R0174S
R0136L	BamHI	50,000 units	R0136L
R0136S	BamHI	10,000 units	NE136S
R3136S	BamHI-HF®	10,000 units	R3136S
R3136L	BamHI-HF®	50,000 units	R3136L
R0143S	BgII	2,000 units	NE143S
R0144S	BgIII	2,000 units	NE144S

R0574S	BsrDI	200 units	R0574S
R0197S	Clal	1,000 units	NE197S
R0176S	Dpnl	1,000 units	NE176S
R0101S	EcoRI	10,000 units	NE101S
R0195S	EcoRV	4,000 units	NE195S
R0104S	HindIII	10,000 units	NE104S
R0124S	HinP1I	2,000 units	R0124S
R0105S	HpaI	500 units	NE105S
R0193S	NcoI	1,000 units	NE193S
R0111S	NdeI	4,000 units	NE111S
R0189S	NotI	500 units	NE189S
R0189L	NotI	2,500 units	NE189L
R0547S	PacI	250 units	NE547S
R0140S	PstI	10,000 units	NE140S
R0157S	SacII	2,000 units	NE157S
R0138S	Sall	2,000 units	NE138S
R0141S	SmaI	2,000 units	NE141S
R0141L	SmaI	10,000 units	NE141L
R0187S	StuI	1,000 units	NE187S
R0145S	XbaI	3,000 units	NE145S
R0146S	XhoI	5,000 units	NE146S

NEB Representative: Bruno Suttles
bsuttles@neb.com
(415) 830-0140

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com