

Tray	NEB #	Buffers	Size	Quantity Available
16	B0518S	Phusion [®] HF Buffer Pack	6 ml	2
Tray	NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Quantity Available
7	E2621S	NEBuilder [®] HiFi DNA Assembly Master Mix	10 reactions	1
Tray	NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Quantity Available
5	M0289S	Antarctic Phosphatase	1,000 units	2
5	M0303S	DNase I (RNase-free)	1,000 units	1
6	M2200S	Quick Ligation [™] Kit	30 reactions	2
5	M0371S	Shrimp Alkaline Phosphatase (rSAP)	500 units	2
6	M0202L	T4 DNA Ligase	100,000 units	2
6	M0202S	T4 DNA Ligase	20,000 units	2
6	M0202T	T4 DNA Ligase	20,000 units	2
6	M0202M	T4 DNA Ligase	100,000 units	2
5	M0201S	T4 Polynucleotide Kinase	500 units	1
7	M0302S	T7 Endonuclease I	250 units	2
Tray	NEB #	Markers & Ladders	Size	Quantity Available
9	N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	2
9	N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	1
10	N0552S	Quick-Load [®] Purple 1 kb DNA Ladder	125 gel lanes (1.25 ml)	1
10	N0550S	Quick-Load [®] Purple 1 kb Plus DNA Ladder	250 gel lanes (1.25 ml)	2
10	N0551S	Quick-Load [®] Purple 100 bp DNA Ladder	125 gel lanes (1.25 ml)	2
SH 3	T2010S	Monarch [®] Total RNA Miniprep Kit	50 preps	out of stock
SH 6	T3050L	Monarch [®] HMV DNA Extraction Kit for Cells & Blood	50 preps	out of stock
SH 5	T2050L	Monarch [®] RNA Cleanup Kit (500 µg)	100 preps	out of stock
SH 2	T1120S	Monarch [®] Spin DNA Gel Extraction Kit	50 preps	out of stock
SH 2	T1120L	Monarch [®] Spin DNA Gel Extraction Kit	250 preps	out of stock
SH 3	T3010S	Monarch [®] Spin gDNA Extraction Kit	50 preps	out of stock
SH 3	T1130S	Monarch [®] Spin PCR & DNA Cleanup Kit (5 µg)	50 preps	out of stock
SH 4	T1130L	Monarch [®] Spin PCR & DNA Cleanup Kit (5 µg)	250 preps	out of stock
SH 1	T1110L	Monarch [®] Spin Plasmid Miniprep Kit	250 preps	out of stock
SH 1	T1110S	Monarch [®] Spin Plasmid Miniprep Kit	50 preps	out of stock
Tray	NEB #	PCR, Polymerases & Amplification Technologies	Size	Quantity Available
11	N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol	2
11	M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	2
14	E3005S	Luna [®] Universal One-Step RT-qPCR Kit	200 reactions	2
14	E3006S	Luna [®] Universal Probe One-Step RT-qPCR Kit	200 reactions	2
15	M3004S	Luna [®] Universal Probe qPCR Master Mix	200 reactions	2
15	M3003S	Luna [®] Universal qPCR Master Mix	200 reactions	2
13	M0480S	OneTaq [®] DNA Polymerase	200 units	2
13	M0480L	OneTaq [®] DNA Polymerase	1,000 units	1
12	M0509L	OneTaq [®] Quick-Load [®] DNA Polymerase	500 units	3

13	M0530L	Phusion [®] High-Fidelity DNA Polymerase	500 units	1
8	E6560S	ProtoScript [®] II First Strand cDNA Synthesis Kit	30 reactions	2
12	M0491S	Q5 [®] High-Fidelity DNA Polymerase	100 units	3
12	M0491L	Q5 [®] High-Fidelity DNA Polymerase	500 units	2
8	M0314S	RNase Inhibitor, Murine	3,000 units	2
11	M0203S	T4 DNA Polymerase	150 units	1
Tray	NEB #	Protein Tools	Size	Quantity Available
10	P7718S	Blue Prestained Protein Standard, Broad Range (11-250 kDa)	150 gel lanes	2
16	P8107S	Proteinase K, Molecular Biology Grade	2 ml	2
Tray	NEB #	Restriction Endonucleases	Size	Quantity Available
1	R3552S	Agel-HF [®]	300 units	2
1	R0558S	Ascl	500 units	2
1	R3136S	BamHI-HF [®]	10,000 units	1
1	R0176S	Dpnl	1,000 units	2
2	R3101S	EcoRI-HF [®]	10,000 units	2
2	R3195S	EcoRV-HF [®]	4,000 units	2
2	R3142S	KpnI-HF [®]	4,000 units	2
2	R3193S	NcoI-HF [®]	1,000 units	2
2	R0111S	NdeI	4,000 units	1
3	R3189S	NotI-HF [®]	500 units	2
3	R3189L	NotI-HF [®]	2,500 units	2
3	R3156S	SacI-HF [®]	2,000 units	2
3	R3138S	Sall-HF [®]	2,000 units	2
4	R3122S	Scal-HF [®]	1,000 units	2
4	R0145S	XbaI	3,000 units	2
4	R0146S	XhoI	5,000 units	1
4	R0180S	XmaI	500 units	2

NEB Representative: Shanna Sawatzki
shanna@neb.com
(206) 889-0401

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com