

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0220S	AluI Methyltransferase	100 units	M0220S
M0298L	Cre Recombinase	250 units	M0298L
M0202T	T4 DNA Ligase	20,000 units	M0202T
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	N3200S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
M0259S	Deep Vent® (exo-) DNA Polymerase	200 units	M0259S
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	M0210S
M0480S	OneTaq® DNA Polymerase	200 units	M0480S
M0481S	OneTaq® Hot Start DNA Polymerase	200 units	M0481S
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	M0530S
M0203S	T4 DNA Polymerase	150 units	M0203S
M0267S	Taq DNA Polymerase with ThermoPol® Buffer	400 units	M0267S
M0254S	Vent® DNA Polymerase	200 units	M0254S
NEB #	Restriction Endonucleases	Size	Stockroom #
R0161S	AccI	1,000 units	R0161S
R3552S	AgeI-HF®	300 units	R3552S
R0507S	ApaI	2,500 units	R0507S
R0136S	BamHI	10,000 units	R0136S
R3136S	BamHI-HF®	10,000 units	R3136S
R0118S	BanI	5,000 units	R0118S
R0144S	BglII	2,000 units	R0144S
R0531S	BsaAI	500 units	R0531S
R0553S	BsiWI	300 units	R0553S
R0502S	BspMI	100 units	R0502S
R0113S	BstXI	1,000 units	R0113S
R0197S	Clal	1,000 units	R0197S
R0176S	DpnI	1,000 units	R0176S
R0101S	EcoRI	10,000 units	R0101S

R3101S	EcoRI-HF [®]	10,000 units	R3101S
R0195S	EcoRV	4,000 units	R0195S
R0139S	HhaI	2,000 units	R0139S
R0104S	HindIII	10,000 units	R0104S
R3104S	HindIII-HF [®]	10,000 units	R3104S
R0155S	Hinfl	5,000 units	R0155S
R0105S	HpaI	500 units	R0105S
R3142S	KpnI-HF [®]	4,000 units	R3142S
R0193S	NcoI	1,000 units	R0193S
R3193S	NcoI-HF [®]	1,000 units	R3193S
R0564S	NgoMIV	1,000 units	R0564S
R0189S	NotI	500 units	R0189S
R0189L	NotI	2,500 units	R0189L
R3189S	NotI-HF [®]	500 units	R3189S
R3189L	NotI-HF [®]	2,500 units	R3189L
R0127S	NsiI	1,000 units	R0127S
R0547L	PacI	1,250 units	R0547L
R0547S	PacI	250 units	R0547S
R0151S	PvuII	5,000 units	R0151S
R3151S	PvuII-HF [®]	5,000 units	R3151S
R3156S	SacI-HF [®]	2,000 units	R3156S
R0157S	SacII	2,000 units	R0157S
R0138S	SalI	2,000 units	R0138S
R3138S	SalI-HF [®]	2,000 units	R3138S
R0569S	SapI	250 units	R0569S
R0187S	StuI	1,000 units	R0187S
R0145L	XbaI	15,000 units	R0145L
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0180S	XmaI	500 units	R0180S
NEB #	RNA Synthesis and Modification	Size	Stockroom #
M0204S	T4 RNA Ligase 1 (ssRNA Ligase)	1,000 units	M0204S

NEB Representative: Allan Waite
awaite@neb.com
(978) 471-0698

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com