

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0289S	Antarctic Phosphatase	1,000 units	M0289S
M0303S	DNase I (RNase-free)	1,000 units	M0303S
M0206S	Exonuclease III (<i>E. coli</i>)	5,000 units	M0206S
M0212S	Klenow Fragment (3'→5' exo-)	200 units	M0212S
M0202S	T4 DNA Ligase	20,000 units	M0202S
NEB #	Glycobiology	Size	Stockroom #
P0705S	PNGase F (Glycerol-free)	15,000 units	P0705S
NEB #	Markers & Ladders	Size	Stockroom #
N3232L	1 kb DNA Ladder	1,000 gel lanes (1 ml)	N3232L
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3200L	1 kb Plus DNA Ladder	1,000 gel lanes (0.5 ml)	N3200L
N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	N3200S
N3231L	100 bp DNA Ladder	500 gel lanes (0.5 ml)	N3231L
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 μmol	N0447S
M0480S	OneTaq® DNA Polymerase	200 units	M0480S
M0481S	OneTaq® Hot Start DNA Polymerase	200 units	M0481S
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	M0530S
M0203L	T4 DNA Polymerase	750 units	M0203L
M0203S	T4 DNA Polymerase	150 units	M0203S
M0273X	Taq DNA Polymerase with Standard Taq Buffer	4,000 units	M0273X
M0273S	Taq DNA Polymerase with Standard Taq Buffer	400 units	M0273S
M0267X	Taq DNA Polymerase with ThermoPol® Buffer	4,000 units	M0267X
NEB #	Protein Tools	Size	Stockroom #
P8107S	Proteinase K, Molecular Biology Grade	2 ml	P8107S
NEB #	Restriction Endonucleases	Size	Stockroom #
R3552S	Agel-HF®	300 units	R3552S
R3136S	BamHI-HF®	10,000 units	R3136S
R0144S	BglII	2,000 units	R0144S
R0540S	BspEI	1,000 units	R0540S
R0113S	BstXI	1,000 units	R0113S

R0176S	Dpnl	1,000 units	R0176S
R3505S	EagI-HF [®]	500 units	R3505S
R3101S	EcoRI-HF [®]	10,000 units	R3101S
R3101T	EcoRI-HF [®]	10,000 units	R3101T
R0195S	EcoRV	4,000 units	R0195S
R3195S	EcoRV-HF [®]	4,000 units	R3195S
R0104S	HindIII	10,000 units	R0104S
R0105S	HpaI	500 units	R0105S
R3142S	KpnI-HF [®]	4,000 units	R3142S
R0525S	MseI	500 units	R0525S
R0193S	NcoI	1,000 units	R0193S
R3193S	NcoI-HF [®]	1,000 units	R3193S
R0111S	NdeI	4,000 units	R0111S
R3131S	NheI-HF [®]	1,000 units	R3131S
R3189S	NotI-HF [®]	500 units	R3189S
R3140S	PstI-HF [®]	10,000 units	R3140S
R0157S	SacII	2,000 units	R0157S
R3138T	SalI-HF [®]	2,000 units	R3138T
R3138S	SalI-HF [®]	2,000 units	R3138S
R0123S	SfiI	3,000 units	R0123S
R0141S	SmaI	2,000 units	R0141S
R3182S	SphI-HF [®]	500 units	R3182S
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0180S	XmaI	500 units	R0180S

NEB Representative: Allan Waite
awaite@neb.com
(978) 471-0698

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com